

OrthoSta

Fixation System

The versatile knotless, crimpless system for stabilizing injured joints.

- ✓ Flat-woven UHMWPE fiber suture tape for high strength and low elongation
- ✓ Novel titanium alloy interference screw for secure, knotless fixation
- ✓ Healix™ antimicrobial technology for infection prevention

PATIENT WEIGHT RECOMMENDATION:

< 40 Lbs (18 Kg) — OrthoSta System

Item#	Description
OSFTL 2036	<i>OrthoSta UHMWPE Woven Tape 2mm, Blunt Needle, 36cm, TI SB-AS, Healix Antimicrobial</i>
OSIS 200.08	OrthoSta Interference Screw, for 2mm Tape, 8mm, T6 Star Recess, TI-6AL-4V
OSTAP 2.0	OrthoSta Tap, for 2mm Tape, Stainless Steel
500035	Screwdriver Blade, Star Driver QCK Shaft, T6
CAN 2.5 160 QCK	2.5mm Cannulated Quick Coupling Drill Bit, 160mm

PATIENT WEIGHT RECOMMENDATION:

> 40 Lbs (18 Kg) — OrthoSta System

Item#	Description
OSFTL 4060	<i>OrthoSta UHMWPE Woven Tape 4mm, Blunt Needle, 60cm, TI SB-PS, Healix Antimicrobial</i>
OSIS 400.13	OrthoSta Interference Screw, for 4mm Tape, 13mm, T8 Star Recess, TI-6AL-4V
OSTAP 4.0	OrthoSta Tap, for 4mm Tape, Stainless Steel
500036	Screwdriver Blade, Star Driver QCK Shaft, T8
CAN 3.2 170 QCK	3.2mm Cannulated Quick Coupling Drill Bit, 170mm



Upcoming Stifle Workshops

Movora Education is proud to offer a comprehensive catalog of hands-on labs and workshops to veterinary professionals worldwide.

Feb.	February 22-23 – VOI Stifle Stabilization: OrthoSta CrCL Technique <i>University of Calgary – Calgary, AB</i>
April	April 13-14 – VOI Stifle Stabilization: OrthoSta CrCL Technique <i>Movora Education Center – Boston, MA</i>
May	May 23-24 – VOI Stifle Stabilization: OrthoSta CrCL Technique <i>Viticus Group, Oquendo Center – Las Vegas, NV</i>
June	June 7-8 – VOI Stifle Stabilization: OrthoSta CrCL Technique <i>Movora Education Center – Boston, MA</i>
July	July 21-22 – VOI Stifle Stabilization: OrthoSta CrCL Technique <i>Viticus Group, Oquendo Center – Las Vegas, NV</i>
Aug.	August 2-3 – VOI Stifle Stabilization: OrthoSta CrCL Technique <i>Movora Education Center – Boston, MA</i>
Dec.	December 12-13 – VOI Stifle Stabilization: OrthoSta CrCL Technique <i>Movora Education Center – Boston, MA</i>

Register to grow your
stifle surgery knowledge
and your practice!

LEARN TODAY!




Learn to Use the OrthoSta Fixation System

Visit education.movora.com for upcoming workshops featuring the OrthoSta. Register to start learning the system today!



Scan Here

To register for an
OrthoSta Lab Today



© 2024 Movora.